

## Professional Pulsing ElectroMagnetic Field (PEMF) Therapy System Curatron 2000XP

- ✓ Pain relief
- ✓ Arthritis
- ✓ Sports Injuries
- ✓ Osteoporosis
- ✓ Bone non-union
- ✓ Backache
- ✓ Carpal tunnel syndrome
- ✓ Neuralgia
- ✓ Epicondylitis
- ✓ Fibromyalgia
- ✓ Fracture healing
- ✓ Pseudo arthrosis
- ✓ Rheumatic pain
- ✓ Tendinitis
- ✓ Low back pain
- ✓ Wound healing



The Curatron 2000XP is the most popular Curatron model for health care professionals and discerning individuals. It is based on state-of-the-art technology generating Frequency Modulated Pulsing ElectroMagnetic Fields (PEMF). These fields are completely safe and give no side effects.

Improved blood supply increases the oxygen pressure, activating and regenerating cells. Optimal peripheral oxygenation is obtained by this professional, high-energy system, resulting in maximum therapy efficacy. Improved calcium transport increases absorption of calcium in bones and improves the quality of cartilage in joints, decreasing pain dramatically. At the same time substantial increase in bone density is obtained, which in turn decreases fracture risk. In addition, bone growth is stimulated causing fractured bones to heal properly.








This technology has undergone many scientific clinical trials worldwide. Tens of thousands of patients have been treated successfully. Of the treated patients, a very high success percentage is reported, creating a dramatic improvement in quality of life. The Curatron 2000XP is the only economical priced computerized therapy device, which uses frequency modulated pulsed electromagnetic field technology and generates sufficient energy to obtain a **real** therapeutic effect.

10 Different treatment programs are available, allowing for the best possible individual automatic parameter settings for the diseases as indicated above. The pulsing electromagnetic fields are applied to the body by means of multiple flat coils, which are located inside the therapy pad. The therapy pad is placed at the area to be treated and is hooked up to the Curatron unit via a cable. The controls of the system are very user-friendly and simple to operate. Select a program number and press the Start button and the complete therapy is performed automatically, under full control of the built-in computer. After the therapy time is over, the energy switches itself off and the unit goes into "stand-by" mode. Treatment can be undergone safely while relaxing or sleeping.

No special training for use of the system is required.  
In order to obtain optimal benefits of the Curatron 2000XP it is to be applied on a regular basis.

**Options:**

-  **Therapy pad of 20" x 28" (50 x 70 cm) emitting up to 30 milli Tesla (300 Gauss)**
  
-  **High Energy Coil: 8" x 8" (20 x 20 cm) for local treatment (knee, elbow, hand, foot) emitting up to 45 milli Tesla (450 Gauss).**
  
-  **Very High Energy Coil: 8" x 8" (20 x 20 cm) for local treatment (knee, elbow, hand, foot) emitting up to 65 milli Tesla (650 Gauss). This option is only available to qualified therapists.**
  
-  **Full body mattress: 26" x 71" (65 x 180 cm) emitting up to 10 milli Tesla (100 Gauss). This mattress is specially designed for full body osteoporosis therapy.**
  
-  **Switch box for convenient selection between 2 different applicators.**